3D Bond™

Graft Binder Cement

Instructions For Use

Augma™
BIOMATERIALS

Table of Contents

(EN)	English	3-10
(FR)	Français	11-18
(RU)	Русский	19-26
(RO)	Romanian	27-34
(IT)	Italiano	35-42
(GR)	Deutsch	43-50
(SP)	Español	51-58
(AR)	العربية	59-66
(HE)	עברית	67-74

Attention: See instructions for use

Do not re-sterilize

Do Not Use if Package is Damaged

Do not re-use

REF

Catalogue number

LOT

Batch code

STERILE R

Sterilized using gamma irradiation

Data of manufacture
(including sterilization)

Use by date

Temperature limitation

Caution:
US federal law restricts this product
for sale by or on the order of a
dentist or physician.

0344

EC REP

MedNet GmbH,
Borkstrasse 10,
48163 Münster, Germany

Augma™
BIOMATERIALS

Manufacturer:
Augma Biomaterials Ltd.
Alon Hatavor 20 St.,
P.O.Box 3089
Caesarea Southern Industrial
Park 3088900, Israel
www.augmabio.com

INSTRUCTIONS FOR USE

3D Bond™ is a Calcium Sulfate graft binder cement. Please read this entire circular before performing the procedure.

This device is for sale by, or on the order of a physician, or licensed practitioner.

COMPOSITION

Each driver contains medical grade Calcium Sulfate.

PRODUCT DESCRIPTION

3D Bond™ is a synthetic, osteoconductive, bioresorbable graft binder cement composed of biphasic Calcium Sulfate in a granulated powder form. 3D Bond™ functions as a scaffold for bone regeneration in dental procedures, and is intended to fill, augment or reconstruct bony defects in the oral and maxillofacial region.

STERILIZATION

3D Bond™ sterilization is carried out using γ -irradiation.

INDICATIONS FOR USE

3D Bond™ is indicated for use in the following ways: by itself in bone regenerative techniques, mixed with other suitable bone filling agents to prevent particle migration in osseous defects, and to provide a resorbable barrier over other bone graft material.

CONTRAINDICATIONS

- The customary contraindications in oral and maxillofacial surgery with other implant materials should be observed.
- Lack of adequate training of the practitioner is a major risk for the success of the implant procedure.

ABSOLUTE CONTRAINDICATION

- Acute and chronic active infection at the site of the implant.

RELATIVE CONTRAINDICATION

- Serious disturbance of bone metabolism
- Serious bone diseases of endocrine etiology
- Severe or difficult to control diabetes mellitus
- Immunosuppressive and radiation therapy
- Ongoing treatment with glucocorticoids mineralocorticoids and with agents affecting calcium metabolism
- Malignancies
- Lactating and pregnancy
- Effect on pediatric patients is not known

PRECAUTIONS AND WARNINGS

- The material is supplied for one time use. DO NOT RE-STERILIZE in case of partially opened or defective primary packing (driver) or secondary packing (peel-off blister). Under these circumstances 3D Bond™ must no longer be used since the sterility of the material is no longer ensured.
- Do not use when the temperature of the product or the saline is below 10°C (50°F). If used in low temperatures, wait until the product restores room temperature (low temperatures will slow the setting reaction of the material).
- It is not recommended to mix 3D Bond™ with blood.
- In large defects, 3D Bond™ can be mixed with granular augmentation materials in a volumetric ratio of 2:1 (for example, mix 1cc 3D Bond™ with 0.5cc other granular material). Other mixing ratios could compromise the graft hardening in-situ. The ratio calculation should be based on the same units (volume or weight).
- The expiration date is printed on the peel-off blister and on the package. Do not use after indicated expiry date (📅).
- 3D Bond™ does not possess sufficient mechanical strength to support load bearing defects prior to tissue in-growth. In cases where load support is required, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes.
- 3D Bond™ must not be used to stabilize screw placement.
- 3D Bond™ is intended for use by clinicians familiar with bone grafting procedures.
- Possible complications are the same as to be expected of autogenous bone grafting procedures. These may include: superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and/or dislodgement, and general complications that may arise from anesthesia and/or surgery.
- Complications specific to oral/dental use are those as may be typically observed for similar bone grafting procedures and may include: tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation.

STORAGE

Store at temperatures between 5°C (41°F) to 30°C (86°F). Avoid contact with a source of heat. Do not store the product in direct sunlight.

HANDLING

3D Bond™ is available as granulated powder packed within a driver. A sterile standard saline (0.9% Sodium Chloride for injection) is required on-site as a setting agent.

IMPORTANT

- Read all steps of the instructions (1-7) before using 3D Bond™.
- We highly recommend practicing the use of 3D Bond™ before first usage.
- At the beginning of the procedure:
 - Reflect the mucoperiosteal flap.
 - Remove the undesired soft tissue from the exposed bone surface.
 - Prepare the defected area for augmentation procedure.

Note that a membrane is not required in most surgical procedures. In addition, primary flap closure is recommended but not essential for proper healing (since soft tissue may grow on top of 3D Bond™).

3D Bond™ can be used also as a membrane above other augmentation materials. In that case, follow the same steps and apply 3D Bond™ over the selected augmentation material in a thickness which is not less than 1mm.

3D Bond™ can be used by mixing with other suitable bone filling agents to prevent particle migration in an osseous defect. Mixing has to be done in a ratio of 2:1 (For example, mix 1cc 3D Bond™ with 0.5cc other granular material). The preferred granule size should not exceed 1.5mm/0.06”.

Optional method (the purpose of the optional method is to reduce bleeding):

- If there is an accumulation or abundant of bleeding at the surgical site, apply 3D Bond™ in two layers.
- The first layer should be thin. Inject a small amount of 3D Bond™ into the defect area and compress it for few seconds with a dry sponge gauze (a significant decrease of bleeding at the surgical site is expected).
- Apply a second and final 3D Bond™ layer above the first one by compressing it using a dry gauze for 10 seconds (the second layer is applied in order to fill the defect cavity entirely).
- You may carve the material as required using a suitable instrument moistened with saline.
- Proceed with soft tissue coverage and wound closure.

STEP BY STEP


Using 3D Bond™ as a Composite Graft

Read all instruction steps (1-7) before using 3D Bond™. It is recommended that the time elapsed from the introduction of saline into 3D Bond™ should not exceed 3 minutes; otherwise the material will lose its pliability.

PREPARATION STEP

Before applying 3D Bond™, make sure you have the items described:

- A** - A sterile 3D Bond™ driver
- B** - Sterile dry gauze pads
- C** - Sterile medical syringe Sterile saline (Sodium Chloride 0.9% for injection)
- D** - Sterile medical needle 21-27G (1.25"-1.5"/ 31.75mm-38.1mm)
- E** - Sterile granules of your preferable bone graft (the preferred granules size should not exceed 1.5mm/0.06")
- F** - Sterile dish
- G** - Sterile spatula

STEP 1.

Add 3D Bond™ to the granules within the dish. Use a volume ratio of 2 parts 3D Bond™ to 1 part granules. (For example, add 1cc 3D Bond™ to 0.5cc granules).

STEP 2.

Use a spatula and mix the dry blend thoroughly.

STEP 3.

Add 0.8 cc saline for one portion of 1.5cc (granules/ 3D Bond™) for complete wettings of the blend (This is important in order to achieve a complete moistening of the blend). (A)

Use a spatula to mix the blend thoroughly until obtaining a homogeneous creamy texture (The texture is ready when it's well mixed and has a shiny and glossy appearance). (B)

STEP 4.


Absorb water surpluses by applying a dry gauze pad to the mixture. Do not use pressure (The texture is ready when it has a matt appearance).

STEP 5.

Use the empty 3D Bond™ driver as a carrier for the putty by pushing the piston to line 1 marked on the driver and re-load it with the putty.



Note: Make sure that the time elapsed from the introduction of saline into 3D Bond™ until it is use, will not exceed 3 minutes.

STEP 6.

Apply the 3D Bond™ composite putty into the required site by pushing the piston and ejecting the putty. (A)

The composite paste should be in direct contact with the bone and slightly overfilled. (B)


STEP 7.

Shape as necessary to achieve the desired form. Compress with a dry gauze pad for 3-5 seconds for condensing the putty.

Remove the gauze pad.

Proceed with soft tissue coverage and wound closure.

STEP BY STEP


Using 3D Bond™ by Itself

Read all instruction steps (1-7) before using 3D Bond™. It is recommended that the time elapsed from the introduction of saline into 3D Bond™ should not exceed 3 minutes; otherwise the material will lose its pliability.

PREPARATION STEP

Before applying 3D Bond™, make sure you have the items described:

- Ⓐ - A sterile 3D Bond™ driver
- Ⓑ - Sterile saline (Sodium Chloride 0.9% for injection)
- Ⓒ - Sterile dry gauze pads
- Ⓓ - Sterile medical syringe
- Ⓔ - Sterile medical needle 21G-27G (1.25"-1.5"/ 31.75mm-38.1mm)

Please note: Lines 1&2 marked on the driver do not indicate quantities or measurements; they represent the handling steps.

STEP 1.


Push the piston to line 1 marked on the driver. Twist and pull to remove the driver's cap.

STEP 2.


Insert the syringe's needle through the driver's head, and eject the saline from the syringe into the 3D Bond™ driver.

- Eject until excess liquid drops out.
- Do not block the driver's aperture

STEP 3.

Hold the driver's head handles firmly and push the piston towards line 2 to expel excess liquid.


STEP 4.

Hold the driver horizontally and remove the driver's head by twisting and pulling.

STEP 5.



Dab the driver aperture on dry gauze for 5 seconds. This will absorb liquid surplus to obtain the desired viscosity.


STEP 6.

Apply the 3D Bond™ paste into the required site by pushing the piston and ejecting the putty. (A)

The paste should be in direct contact with the bone and slightly overfilled. (B)

**STEP 7.**

Shape as necessary to achieve the desired form. Apply a dry gauze pad for 3-5 seconds and condense the putty. Remove the gauze pad. Proceed with soft tissue coverage and wound closure.

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